



Kimal is proud to announce a new addition to the SurgiGUARD<sup>®</sup> family. SurgiGUARD<sup>®</sup> G is a range of sterile packed gowns that offers quality and protection to meet your clinical needs.

All gowns conform fully and meet the standard performance for critical product area requirements of the EN13795 standard.



## A range of non reinforced SSMMS gowns offering moderate protection

- SSMMS is fluid repellent
- Double layer of meltblown material improves performance
- Packed with two sterile hand towels
- All gowns feature an overlap sterile back with adjustable velcro fastening
- Soft fabric cuff for extra comfort

Description	Code	Unit of Issue
Standard SSMMS gowns with velcro fastening - Large	K62/G01	30
Standard SSMMS gowns with velcro fastening - Extra large	K62/G02	30

## A range of non reinforced Spunlace gowns offering moderate protection

- Spunlace is fluid repellent
- Packed with two sterile hand towels
- All gowns feature an overlap sterile back with adjustable velcro fastening
- Soft fabric cuff for extra comfort

Description	Code	Unit of Issue
Standard spunlace gowns with velcro fastening - Large	K62/G03	30
Standard spunlace gowns with velcro fastening - Extra large	K62/G04	30

